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Prestomat[®]

compact

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

Umschlags. 2

Vakat !

Prestomat® compact

D	Gebrauchsanweisung	5
GB	Instructions for use	23
F	Mode d'emploi	39
I	Istruzioni per l'uso	55
E	Instrucciones para el uso	73

Seite 4

Vakat !

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Seite 22

Vakat !

1	Content	
2	Proper use	24
2.1	Scope of delivery	24
2.2	Available accessories	24
2.3	Technical data	24
2.4	Safety information	25
3	Putting into operation	25
3.1	Compressed air connection	25
3.2	Connection and assembly of the Prestomat Printer	26
4	Description of the individual components	26
4.1	Operation panel	27
4.2	Displays	28
4.3	Display of casting sequence	28
4.4	Control keys	29
4.5	Functions and adjustment ranges	30
5	Programming	31
5.1	Programming of the casting sequence	31
5.2	Programming the increase of the melting temperature	32
5.3	Program tables	32
6	Working instructions	33
6.1	Selection of the program	33
6.2	Melting and casting	33
7	Errors	35
8	Care and maintenance	36
8.1	Control of the casting chamber rotation and the vacuum	36
8.2	Values of the simulation programs	36
8.3	Cleaning of the furnace chamber	37
9	Statement of conformity	38

2 Proper use

The Prestomat® compact is an electrically driven vacuum pressure casting unit for melting and casting of dental alloys. A newly developed ceramic thermocouple allows a max. heating temperature of 1500 °C. The motor-driven, automatic turning mechanism of the casting chamber contributes to achieve reproducible casting results at any time. Numerous functions that can be easily controlled ensure operational convenience. Safe operating is achieved through a control which excludes malfunctions. By connecting the Prestomat Printer, casting processes can be logged.

Carefully read these instructions for use prior to the installation and initial operation. The information in these instructions for use must be observed since operation that is not in compliance with these instructions may include hazards for the operator.

2.1 Scope of delivery

- 1 unit complete with vacuum pump and compressed air connection
- 1 eccentric ring
- 1 crucible holder
- 1 crucible stand
- 1 cleaning insert
- 1 tweezers
- 1 instructions for use/statement of conformity
- 1 guarantee card
- 1 pressure tank certificate

REF 5522 0001

2.2 Available accessories

- Prestomat Printer REF 5522 0200
- graphite crucible Pc (package with 10 pcs)
REF 5522 0101
- ceramic crucible Pc (package with 10 pcs)
REF 5522 0102

2.3 Technical data

Voltage supply:	230 V~, 50/60 Hz
Maximum power consumption:	1,5 kW
Dimensions:	300 mm (W) x 420 mm (H) x 345 mm (D)
Weight:	approx. 390 N (40 kg)
Max. heating temperature:	1500 °C
Compressed air connection:	0,5–1,0 MPa
Casting programs:	16 (0 – 9, A – F)

The unit is designed for indoor use only. The room temperature should be between 5 °C and 40 °C. Main supply voltage variations not larger than $\pm 10\%$ of the nominal value. Contamination degree 1 or 2 according to IEC 664 / section 3.7.3.

2.4 Safety information

Prestomat compact must only be operated with original spare parts and accessories. Only in this way the performance data provided can be achieved and the required operational safety is ensured.

Hazardous voltages are fed inside the unit. Maintenance work and repairs must only be carried out by service technicians authorized by DeguDent.

The unit must be protected against moisture in order to avoid short circuits and damage to the unit. Do not expose the unit to water vapor.

The chamber is filled with a pressure of approx. 0.3 MPa. The casting chamber must not be opened while it is under pressure. Do not touch the casting chamber during the casting process.

Do not look inside the melting furnace if the temperature is above 800 °C. For this purpose the viewing glass that is included in the scope of delivery must be used.

Unplug the unit if it is not to be used for longer periods.

3 Putting into operation

Place the unit on a flat, solid base with a minimum distance of 10 cm to the wall or other objects. Do not place the unit on flammable surfaces. Do not place any flammable objects near the unit. The crucible stand that is included is suitable to store hot crucibles. Keep the unit away from dusty locations.


Connect the unit to the main supply. Be sure that the power supply corresponds to the electrical specifications of the unit.

3.1 Compressed air connection

At the rear of the unit is a pressure reducer with a compressed air inlet (B) for nylon or polyurethane hoses with an outside diameter of 8 mm. Connect the unit to the compressed air supply with a pressure of 0.5–1.0 MPa. The hose connection should be short. In order to avoid hazards only compressed air hoses that are suitable for the working pressure must be used. The unit features a rapid action hose coupling. The compressed air hose must be plugged into the rapid action hose coupling to the stop. Check the correct position by quickly pulling back the hose.

If moisture has been accumulated in the separator (C), press the discharge valve.


The pressure in the casting chamber (H) has been set to 0.3 MPa (± 0.03 MPa). Pressure that is too low results in defective castings. Pressure that is too high may cause leakage. At a pressure of more than 0.4 MPa the safety valve (A) is released.

- A) Safety valve
- B) Compressed air connection
- C) Moisture separator

3.2 Connection and assembly of the Prestomat Printer


The serial port (P) for the connection of the printer is located at the right side of the housing. The printer is attached into the openings on the right side of the housing and fixed with the screwing located at the printer. For programming of the printer please refer to the Prestomat Printer operating instructions.


- ! Caution:**
Power is supplied to the printer via the serial port. Therefore do not connect any commercially available printer to the Prestomat compact since this may result in damage to the units !

4 Description of the individual components

In order to open the lid of the casting chamber (F), move the locking lever upward. In order to close the lid, move the locking lever (G) downward until it locks.


- D) Operation panel
- E) Pressure display
- F) Lid of casting chamber
- G) Locking lever
- H) Casting chamber
- I) On/off switch

- J) Eccentric ring
- K) Crucible holder
- L) Crucible

4.1 Operation panel

4.2 Displays

4.2.1 Temperature Display (M)

Displays the actual temperature of the casting chamber. If required, it can also display the stored information of the casting programs as well as error messages. With the help of the corresponding control lamp that is lit it is possible to recognize whether the value that is displayed is a temperature (°C) or time (⌚) value.

4.2.2 Program display (N)

Shows the corresponding program out of a total number of 16 possible programs (0 – 9, A – F) that is started or running.

4.2.3 Preheating

In order to obtain a uniform temperature inside the casting chamber and thus optimum conditions for the following casting processes, the unit automatically heats up to 1300 °C when it is switched on. The control lamp flashes during the preheating process. When the preheating temperature is reached, the temperature decreases to the starting temperature that has been entered and the control lamp is lit continuously. If a melting temperature of 0 °C has been entered, the unit recognizes no entry and will not start the preheating program either.

4.3 Display of casting sequence (O)

Shows the sequence of the casting process. Flashing control lights indicate which process is just being run; control lamps that are continuously lit indicate the completion of the corresponding process. Control lamps that are not lit indicate that this process has not been carried out or that it has already been completed.

4.3.1 Starting temperature (Start °C)

The display flashes as long as the starting temperature has not been reached. If the starting temperature is reached, the control lamp is lit continuously. It is also lit during programming and confirming of the melting temperature. After finishing the melting process and the casting chamber is open, the control lamp is flashing again until the starting temperature is reached.

4.3.2 Melting temperature

Until the melting temperature is reached, the control light is flashing. If the melting temperature is reached, the control lamp is lit continuously. It is also lit during programming and confirming the melting temperature. It goes out when the casting process is finished and the lid of the casting chamber is open.

4.3.3 Time of further heating (Further heating sec.)

The time of further heating is started by pressing a button (Gießen). The control lamp flashes until the end of this time. After the end of the time of further heating it is lit continuously and an acoustic signal can be

heard. If the casting chamber is closed, it will start to turn into the casting position after an evacuation time of approx. 20 seconds.

If a time of further heating of 15 min is exceeded, "FAIL" is shown in the display. The current supply to the thermocouple is automatically disconnected. In order to reactivate the current supply, press any key on the operation panel. Then the temperature will level out at the starting temperature that has been programmed.

4.3.4 Pressure delay time (Pressure delay sec.)

Depending on the type and the quantity of alloy a pressure delay time between 0.0 and 0.3 seconds can be preselected. In case of a higher viscosity of the melted alloy an extended pressure delay time is required. If a large alloy quantity is melted, it is recommended to select an extended pressure delay.

During the casting process, i.e. if the casting chamber is turned downward by 180° and the vacuum pump is working, the control lamp is lit. It is lit when the casting process is completed and during programming and confirming of the casting time. If the casting process is completed and the casting chamber is opened, it goes out. Enter the required time until the alloy solidifies.

4.3.5 Cooling time

After the casting process is finished the control lamp is flashing. It lit's continuously while the casting chamber moves back to the starting position and during programming and confirming the cooling time. After the casting process is finished and the lid of the casting chamber is open it goes out.

4.4 Control keys

Press this key in order to enter into the program mode. By pressing this key again you will leave this mode. Programming is only possible during stand-by operation and during the melting process. During the casting process only the confirmation is possible.

Press this key to change to a new program. If it is pressed once, you will move up (▲) or down (▼) one program. Press this key in the programming mode to change the values in the individual programs.

Press this key to store the values that have been entered and to move to the next program. The programmed data are confirmed with this key. It only works in the programming mode. No values will be changed if this key is pressed.

Press this key to start the melting process. The temperature will rise to the melting temperature. While the program is running, the diode of the melting temperature is on. If this key is pressed again, the melting process is aborted. If a melting temperature of 0 °C has been programmed, the unit is unable to handle the program; program course and the key will not work.

Press this key, if the melting temperature that has been entered is insufficient. The feature of the melting temperature increase allows for increase of the programmed melting temperature step-by-step (upon delivery the unit is programmed at **+10 °C**). Repeated pressing of this key results in a further increase of the melting temperature. This step can be repeated as often as desired until the maximum temperature is reached. The diode flashes during the heating process. It is lit continuously if the additional melting temperature is reached. The melting temperature that has been entered is not increased and remains unchanged during a subsequent casting process.

The desired temperature increase can be entered individually (see item 5.2 “Programming the melting temperature increase”).

If this key is pressed, the casting chamber will vibrate for a short moment. This helps to recognize more clearly whether the alloy has been completely melted or not. This key also works during the melting process.

Do not look inside the casting chamber at temperatures above 800 °C. Use the shaded viewing glass that is included !

After the alloys has been melted completely, the entered time of further heating is activated by pressing this key. After the end of the time of increased heating an acoustic signal can be heard. If the casting chamber is closed, evacuating is initiated. After the corresponding vacuum has been reached, the casting process is started. Please remember that approx. 20 seconds are required for evacuating the casting chamber. This must be considered when programming the time of further heating.

Example: If the manufacturer of the alloy recommends a time of increased heating of 60s, the time of increased heating in the corresponding casting program should be programmed to 40s.

4.5 Functions and adjustment ranges

Functions	Adjustment ranges	Smallest adjustment step
+ °C	1–50 °C	1 °C
Starting °C	0–1300 °C	1 °C
Melting °C	0–1500 °C	1 °C
Time of increased heating sec.	0 s – 7 min 30 s	5 s
Pressure delay sec.	0.0–0.3 s	0.1 s
Cooling	5 s – 10 min	5 s

5 Programming

Programming allows to store the values of 5 functions in the casting sequence. Additionally, an increase of the melting temperature (see item 5.2) can be stored. The programmed melting temperature increase is used by all programs.

If the unit is switched on, it is in the standby-mode and adjusts the starting temperature of the program that has been selected.

All steps of the casting sequence can be programmed individually. Please enter these steps into the table under item 5.3 in order to be able to program them a new if they are lost.

- Starting temperature
- Melting temperature
- Time of further heating
- Pressure delay
- Cooling time

The values depend on the processing instructions of the respective alloy. Therefore casting programs for 16 different alloy types can be preprogrammed. Upon delivery the programs 0–9 have been preprogrammed. The corresponding values have been matched with 10 different DeguDent alloys, which can be seen in the table under item 5.3. If required, these values can also be reprogrammed individually. The program items A–F have not been preprogrammed.

5.1 Programming the casting sequence

Step	Key	Function
1	▲ ▼	Selecting a program (16 programs, 0–9 or A–F).
2	Progr.	Lamp is on during programming. Control lamp of starting temperature is lit and the stored starting temperature is shown in the display.
3	▲ ▼	Select the desired value by pressing ▲ or ▼.
4	E	Stores the entered value and moves to the next program item. The control lamp of the next program item is lit and shows the stored value in the display.

If steps 3 and 4 are repeated, all necessary data of the respective program items will be stored.

Pressing the “Progr.” key deactivates the programming mode and the unit returns to the standby mode.

Programming is not possible during a casting process. By pressing the “Progr.” and “E” keys the programmed values can be fetched during the casting sequence. The casting sequence is not disturbed through this step.

5.2 Programming the melting temperature increase

Step	Key	Function
1	Progr.	Activates the programming mode (access to all programs).
2	+ °C	When pressing this key the lamp is lit and the stored temperature is shown in the display.
3	▲ ▼	Select the desired value by pressing ▲ or ▼ .
4	E or Progr.	Stores the data, ends the programming mode and returns to the old position.

The adjusted value of the melting temperature increase is used by all programs.

5.3 Program tables

Melting temperature increase (upon delivery + 10 °C)

Table 1

Program number (alloy)	Starting °C	Melting °C	Time of further heating	Pressure delay	Cooling
0 Degunorm®	850	1150	15 s	0.1 s	1 min
1 Degunorm® supra	850	1240	45 s	0.1 s	1 min
2 Degunorm® pur	850	1200	45 s	0.1 s	1 min
3 BiOcclus® HT	1000	1280	45 s	0.2 s	1 min
4 Degudent® U	1100	1400	2 min 45 s	0.2 s	1 min
5 Degudent® H	1100	1400	1 min 45 s	0.2 s	1 min
6 Degulor® M	900	1200	15 s	0.1 s	1 min
7 Deva® 4	1150	1450	2 min 45 s	0.2 s	1 min
8 Stabilor® NF IV	950	1180	15 s	0.1 s	1 min
9 Degupal® G	1100	1450	2 min 45 s	0.2 s	1 min 30 s

Table 2

Program (alloy)	Starting °C	Melting °C	Time of further heating min/s	Pressure delay s	Cooling min/s
A					
B					
C					
D					
E					
F					

6 Working instructions

Start the unit with the main switch. “CAST” is shown in the display when the casting chamber is open. When the lid is closed, an acoustic signal can be heard, “DOOR” is shown in the display and the unit does not heat up.

6.1 Selection of the program

When the unit is switched on, the casting program is displayed that was used last. By pressing the keys ▲ or ▼ any program can be selected. The “Vorwärmen (preheating)” control lamp flashes. In order to achieve uniform distribution of heat in the casting chamber, the thermocouple heats up the casting chamber to 1300 °C and then cools down to the adjusted starting temperature (Starting °C). The necessary preheating time which depends on the programmed starting temperature requires approx. 5–10 min.

When the preheating process is completed, the control lamps of the preheating and the starting temperature are lit. An acoustic signal can be heard.

6.2 Melting and casting

Place the crucible recommended for the alloy (graphite or ceramic) into the furnace chamber, fill the alloy into the crucible and lock the crucible securing device.

Only original crucibles must be used for the Prestomat compact

(see item 2.2 Available accessories).

Press the “Schmelzen (melting)” key. Wait until the melting temperature has been reached and the alloy has been melted completely. If the preprogrammed melting temperature is insufficient, it can be increased

by the programmed temperature (upon delivery 10 °C) by pressing the + °C key. The time of increased heating is activated by pressing the "Gießen (casting)" key. The preprogrammed time value starts to count down.

Place the preheated casting ring with the crucible facing downward into the casting chamber shortly before the end of the time of increased heating and move the eccentric ring clockwise until the casting ring is located in the center. If the casting ring is larger than 6 x, remove the eccentric ring and center the ring on the opening of the crucible.

The eccentric ring gets hot during the melting process. The ring must only be touched with tongs !

After the end of the time of increased heating an acoustic signal can be heard. Close the lid of the casting chamber and lock it with the locking lever. Be sure that the lever has been pressed up to the stop.

If the chamber has been locked correctly, it is automatically evacuated. After sufficient evacuation of the casting chamber, an acoustic signal can be heard and the chamber automatically turns by 180°.

In order to avoid injuries, do not touch the casting chamber any longer after the acoustic signal has been heard. Never try to open the casting chamber during the casting process !

If the rotating movement of the casting chamber is blocked, it turns back to its starting position after approx. 5 s. Err 3 or Err 4 is shown in the display (see item 7).

After the end of the preprogrammed pressure delay time the alloy is cast under pressure. After the end of the programmed cooling time the casting chamber turns back to its starting position and an acoustic signal can be heard. Open the lid of the casting chamber by unlocking the locking lever. Turn the eccentric ring anticlockwise using the crucible tongs, remove the casting ring and place it on a heat-resistant surface. Remove the crucible and put it into the crucible stand.

After a completed casting process the casting chamber, eccentric ring, casting ring and crucible are hot. Do not touch because of the risk of getting burned !

The unit switches back to the standby mode. The temperature levels out at the programmed starting temperature in order to be ready for the next casting process. Leave the lid open. If the lid is closed, an acoustic signal can be heard and "DOOR" is shown in the display. The unit will not heat up then.

7 Errors

Message	Error	Causes and remedial measures
Err 2	Processor error	Switch off the unit. If the same error reoccurs when the unit is switched on again, disturbing pulses from the power supply line resp. from other consumers may be the cause. Another circuit should be chosen to solve the problem.
Err 3 Err 4	Incorrect rotation of the casting chamber	Switch off the unit. Check whether the rotation of the casting chamber is blocked by an obstacle. If an obstacle can be excluded, please inform a service technician.
Err 5	Thermocouple error	Switch off the unit. Thermocouple defective. Have the thermocouple exchanged by the service.
Err 6	Thermosensor error	Switch off the unit. Thermosensor defective. Have the thermosensor exchanged by the service.
Err 7	Vacuum or vacuum control system error	Switch off the unit. Chamber pressure does not reach the diminished pressure required for casting. Clean the washer at the chamber lid. If required, inform a service technician.
Err 8	Storage loss	Switch off the unit. If the unit has not been used for a longer period, it will lose the data that have been programmed individually. If the unit is switched on again, the programs are set back to the values adjusted upon delivery. Therefore the individual program values should be entered into the table under item 5.3. Before new programming is carried out, the unit must be left on for approx. 60 min so that the storage battery is recharged.
Err 9	Incorrect heating	Switch off the unit. An incorrect signal during the heating phase has been detected. If Err 9 is no longer displayed when switching on the unit, this may be due to the electric current supply. Please inform a service technician.
Fail	A time of further heating of 15 min was exceeded	Press any key on the operation panel. The temperature will level out at the preprogrammed starting temperature.

8 Care and maintenance

Care and maintenance include works which ensure trouble-free operation of the Prestomat compact.

The unit must be regularly checked for mechanical damage. In this way accidents and unit failures can be prevented.

Contaminations on all movable parts must be removed immediately. Only in this way smooth running and perfect function are ensured.

The unit features an internal power supply via a storage battery so that data that are stored in the computer control will not be lost in case of a power failure. During operation the storage main supply. The storage battery is maintenance-free.

Hazardous voltages are fed inside the unit. Therefore maintenance work and repairs must only be carried out by service technicians authorized by DeguDent.

Manipulations of the unit reduce the operational safety and include the risk of damage to persons and materials.

8.1 Control of the casting chamber rotation and the vacuum

In order to obtain perfect, reproduceable casting results, a regular control of the inverting movement of the casting chamber and the vacuum level is recommended. For this purpose the Prestomat compact features simulation programs under the program numbers 0–2. In the simulation program the thermo-couple is not heated up.

Proceed as follows:

While pressing the “E” key, switch on the unit. By pressing the “Schmelzen (melting)” resp. “Gießen” keys a casting process is simulated. Check the proper function of the programs.

8.2 Cleaning of the furnace chamber

The cleaning insert included in the scope of delivery is suitable for cleaning the furnace chamber. Cleaning should only be performed when the furnace chamber is cold.

Remove the eccentric ring from the furnace chamber. Instead of the ring place the cleaning insert with the v-shaped opening facing upward into the casting chamber. Close the lid of the casting chamber and manually turn the casting chamber into the casting position. In this way the dirt is collected by the cleaning insert. Then turn the casting chamber back to the starting position, open it and remove the cleaning insert. Place the eccentric ring back into the casting chamber.

In order to avoid burns, please remember that the cleaning insert heats up quickly inside the casting chamber and must therefore not be touched !

In order to prevent contamination within the furnace chamber, the lid of the casting chamber should always be closed when the unit is switched off.

8.3 Values of the simulation programs

Program number	Starting °C	Melting °C	Time of further heating min/s	Pressure delay s	Cooling min/s
0	1000	1130	10s	0.3s	10s
1	1100	1435	10s	0.3s	10s
2	850	1000	10s	0.1s	10s

9 Declaration of conformity

EU Declaration of Conformity, EU Machines Directive (98/37/EG)

DeguDent GmbH

Postfach 13 64 · 63403 Hanau
www.degudent.de

We hereby declare that the:


Prestomat® compact / vacuum pressure casting unit

conforms with the basic requirements of the following EU directives:

- 1. EU Machinery Directive 98/37/EC**
- 2. EEC Directive on electrical equipment designed for use within certain voltage limits 73/23/EEC**
- 3. EEC Directive relating to electromagnetic compability 89/336**

The following harmonized standards were employed:

EN 292-1:1991, EN 292-2:1995, EN 60204-1:1998
EN 61010-1:1994, EN 50081-1:1993, EN 50082-1:1997

Dr. Alexander Völcker
Research & Development

Torsten Schwafert
Director business area
dental laboratory products

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